

RAYPLAN 10A

System Environment Guidelines



RayPlan
RayStation

10A

Traceback information:
Workspace Main version a576
Checked in 2020-06-05
Skribenta version 5.3.050

Disclaimer

For information on functionality not available for regulatory reasons, see the Regulatory Information in the RayPlan Instructions for Use.

Declaration of conformity



Complies with 93/42/EEC Medical Device Directive as amended by M1 to M5. A copy of the corresponding Declaration of Conformity is available on request.

Copyright

This document contains proprietary information that is protected by copyright. No part of this document may be photocopied, reproduced or translated to another language without prior written consent of RaySearch Laboratories AB (publ).

All Rights Reserved. © 2020, RaySearch Laboratories AB (publ).

Trademarks

RayStation, RayBiology, RayCare, RayStore, RayCloud, RayPlan, RaySearch Laboratories and the RaySearch Laboratories logotype are trademarks of RaySearch Laboratories AB (publ).

Third-party trademarks as used herein are the property of their respective owners, which are not affiliated with RaySearch Laboratories AB (publ).

RaySearch Laboratories AB (publ) including its subsidiaries is hereafter referred to as RaySearch.

TABLE OF CONTENTS

1	INTRODUCTION	7
1.1	Contact information	8
1.2	Definitions, acronyms and abbreviations	8
1.3	References and related documents	8
2	HARDWARE ENVIRONMENT	9
3	RAYPLAN MACHINES	11
3.1	RayPlan Client Machines	12
3.2	Stand-alone machines	13
4	SQL SERVERS	15
4.1	SQL Server specifications	16
5	IT CONSIDERATIONS AND MAINTENANCE	19
5.1	General	20
5.2	RayPlan storage SCP server	20
5.3	Third party software	20
5.4	Active Directory (AD)	21
5.5	AD Organizational Unit	22
5.6	SQL	22
5.7	Licensing of third party software	22
5.8	RAM	23
6	MICROSOFT REFERENCE INFORMATION	25

1 INTRODUCTION

About this guideline

This document provides a description of the primary requirements of the hardware and additional software required to run RayPlan. A summary of the ongoing maintenance tasks that should be considered from an IT management perspective is also included.

Any recommendations in this guideline shall be in adherence with *RSL-D-RP-10A-IFU, RayPlan 10A Instructions for use*. If there is any doubt regarding a configuration please contact RaySearch Laboratories.

In this chapter

This chapter contains the following sections:

1.1	Contact information	p. 8
1.2	Definitions, acronyms and abbreviations	p. 8
1.3	References and related documents	p. 8

1.1 CONTACT INFORMATION

Please raise any questions with your local RaySearch Laboratories representative or contact RaySearch Support on support@raysearchlabs.com.

1.2 DEFINITIONS, ACRONYMS AND ABBREVIATIONS

Definitions, acronyms and abbreviations	Meaning
AD	Active Directory
CPU	Central Processing Unit
CUDA	Compute Unified Device Architecture. A computing platform and API developed by NVIDIA for Graphics Processing Units.
DR	Disaster Recovery
ECC RAM	Error Correcting Code Random Access Memory
GPU	Graphics Processing Unit
HDD	Hard Disk Drive
LAN	Local Area Network
PACS	Picture Archiving and Communication System
SAN	Storage Area Network
SSD	Solid State Drive
VLAN	A "section" of a Local Area Network exclusively for traffic of a certain type or a certain application.
WAN	Wide Area Network

1.3 REFERENCES AND RELATED DOCUMENTS

The documents referenced in the following list can be found in **Documentation** that can be launched from the RayPlan Launcher or in the RayPlan Deployment Package.

- *RSL-D-RP-10A-IFU, RayPlan 10A Instructions for use*
- *RSL-D-RP-10A-DCS, RayPlan 10A DICOM Conformance Statement*
- *RSL-D-RP-10A-SEAT, RayPlan 10A System Environment Acceptance Test Protocol*

2 HARDWARE ENVIRONMENT

Introduction

RayPlan operates within a Microsoft Windows environment with Microsoft SQL Server hosting information relating to the clinic, machines and patients in separate databases. The RayPlan client machine obtains the relevant information from the databases. There is also an additional resource database which defines RayPlan Systems. A System is made up of a single clinic and machine database and one or more patient databases.

Additionally, information can be queried, imported or otherwise from application entities (e.g., PACS) using the DICOM format (see details within the *RSL-D-RP-10A-DCS, RayPlan 10A DICOM Conformance Statement*). This information is saved to the relevant database(s).

RayPlan Client at desk

The RayPlan software is installed on a machine that is located at the users' "desk".




Figure 1. The RayPlan Client at desk. **Note that 64 GB RAM is preferred where possible. RAM DIMMs should be installed to provide optimal RAM usage, i.e., units of 6 for Intel Xeon Scalable (48/96 GB) or units of 4 for Intel Xeon E5 (32/64 GB).*

3 RAYPLAN MACHINES

Introduction

The sections in this chapter should be considered in relation to physical RayPlan workstations.

Note: *For performance, update and operation reasons we highly recommend that only RayPlan related software is installed on a client machine.*

In this chapter

This chapter contains the following sections:

3.1	RayPlan Client Machines	p. 12
3.2	Stand-alone machines	p. 13

3.1 RAYPLAN CLIENT MACHINES

The specifications in this section apply to machines being used to host the RayPlan Client software.

Processor(s)

For RayPlan systems that will not include RayArc (VMAT) operations a single 6 core, 3.0GHz or higher processor is sufficient. Where RayArc (VMAT) is to be used a 2nd processor of similar specification should be included. The latest CPUs in the Intel Xeon Scalable Processor series should be considered when investing in new hardware.

Graphics card (GPU) for graphics rendering

An NVIDIA graphics card with a minimum of 1GB RAM is required for the basic rendering operations within RayPlan 10A. For enterprise operations, RaySearch recommends the NVIDIA Quadro cards, principally based on reliability and recommendations from NVIDIA.

Graphics card (GPU) for accelerated computations

In certain RayPlan 10A configurations RayPlan can run certain aspects of the software computations on GPU instead of the CPU. The primary advantage with this is performance.

The graphic card/operating system combinations populated with a GPU driver version in the following table have been verified for GPU accelerated computation operations. **Bold** options are recommended for new purchases.

Graphic cards	Windows 10
NVIDIA RTX8000	441.12
NVIDIA RTX6000	441.12
NVIDIA RTX5000	441.12
NVIDIA P6000	441.12
NVIDIA P5000	441.12
NVIDIA M6000 24 GB	441.12
NVIDIA M6000 12 GB	441.12
NVIDIA M5000	441.12

Note: *Additional validations may have been performed. Certificates are available from support@raysearchlabs.com.*

Please see reference information within the RSL-D-RP-10A-IFU, RayPlan 10A Instructions for use regarding GPU implementation.

Display

The recommended resolution for the display used together with RayPlan is 1920 x 1200 (preferred) or 1920 x 1080 pixels.

RAM

32/48 GB of DDR3 ECC RAM. ECC RAM is required for system memory in all machines. RAM DIMMs should be installed to provide optimal RAM usage, i.e., units of 6 for Intel Xeon Scalable (48/96 GB) or units of 4 for Intel Xeon E5 (32/64 GB).

HDD

The main considerations with respect to HDD for client machines are: hosting the OS (including swap file, drivers, etc.), RayPlan client software, and potentially patient exports.

Local drives are recommended for client machines. Where possible RAID 1 (or potentially higher) should be implemented with SSD drives, if practical. 150 GB per drive or higher is recommended.

Network

1 Gb/s Ethernet is required, and where possible 10Gb/s connectivity between SQL Server and clients would be beneficial. Where possible, it is recommended that the link between the RayPlan clients and the SQL server is within a separate VLAN or similarly segregated network to avoid contention issues.

OS

It is recommended to use Windows 10 Professional or higher to host RayPlan 10A client software.

Windows 10 Professional should be a minimum of version 1709.

Note: *The possibility of running RayPlan on Microsoft Server 2019 is intended to be assessed during 2020.*

3.2 STAND-ALONE MACHINES

Where the machine will act as a stand-alone machine, i.e., host the SQL databases and RayPlan, additional processor power, RAM and hard drive space can be considered.

4 SQL SERVERS

Introduction

If the present virtual infrastructure has the necessary resources, capacity and connectivity, then it can be considered for the hosting of the SQL Server instance(s). The performance and connectivity of the SQL Server host is a critical factor in the performance of opening, saving and closing patients.

SQL tools are available to assess the viability of SAN resources for hosting the RayPlan Databases. Traffic will also be generated for Auto Recovery. The Auto Recovery feature can be turned off.

Using the RayPlan Storage Tool, databases can be moved to other database servers and remain visible to end users via the indexing service.

4.1 SQL SERVER SPECIFICATIONS

Software

As a guideline for storage, 1 GB per patient can be taken as a reasonable indicator. This may change with 4DCT image sets and is also directly related to the nature and number of image sets added during treatment.

Processor(s)

Single 4 core 3.47 GHz, Xeon CPU or better.

The key aspects in the specification decision for the processor on the SQL Server are data upgrades and SQL maintenance. As RayStorage can reduce the burden of upgrading patients, the planned approach can impact the processor drive requirements for the SQL Server.

Graphics Card

The graphics card on the SQL Server is not important for the operation of SQL Databases, however a basic NVIDIA Graphics card may be considered for DR purposes.

RAM

32 GB of DDR3 ECC RAM. RAM DIMMs should be installed to provide optimal RAM usage, i.e., units of 6 for Intel Xeon Scalable (48/96 GB) or units of 4 for Intel Xeon E5 (32/64 GB).

RAM requirements depend mainly on the profile of patients.

HDD

Entry point: OS: 200 GB SSD (RAID 1) Data: 2 TB or higher (RAID 10).

The numbers and complexities of patients may require a higher level of SQL Server or increased hard drive space to store data files, log files and file streams. This will need to be reviewed as part of the planning process. If available, SSDs can be considered for the hosting of Operating System and (full or partial) storage at all levels. Resilience is a key factor in the decision making process. 64k Block size formatting is recommended by Microsoft for SQL data storage drives. The SQL temp database should also be hosted on SSD where possible.

Network

Minimum 1 Gb/s Ethernet (With teaming considered).

It may be beneficial to team one or more network cards specifically for the SQL Server Service, multiple network cards would then be required. Please note that network infrastructure specifications may restrict the ability to team NICs.

Where possible, backups and any other non-RayPlan client related activity should not use the same communications routes as RayPlan Client communications.

OS

Windows Server 2016 or Windows Server 2012 R2 can be used for hosting the RayPlan SQL software.

SQL Server Software

RayPlan 10A is validated for Microsoft SQL Server 2016 SP2 hosting the databases. Microsoft SQL Server 2017, Microsoft SQL Server 2014 SP3 and Microsoft SQL Server 2012 SP4 can also be used, but may not support integration installations with other RaySearch products. In general Microsoft only sells the latest version – with downgrade rights, but where there is a choice RaySearch recommends the purchase of the latest version with downgrade rights. Purchase of the latest version is also encouraged as it appears that operational performance is an area of focus for development, for Microsoft, in relation to SQL.

Note: *MSSQL Server 2012 may no longer be supported in upcoming releases of RayPlan.*

5 IT CONSIDERATIONS AND MAINTENANCE

Introduction

This section describes the IT considerations and maintenance for RayPlan and the RayPlan system environment.

Note: *For performance, update and operation reasons we highly recommend that only RayPlan related software is installed on a client machine.*

In this chapter

This chapter contains the following sections:

5.1	General	p. 20
5.2	RayPlan storage SCP server	p. 20
5.3	Third party software	p. 20
5.4	Active Directory (AD)	p. 21
5.5	AD Organizational Unit	p. 22
5.6	SQL	p. 22
5.7	Licensing of third party software	p. 22
5.8	RAM	p. 23

5.1 GENERAL

Power options

The Microsoft Power Options (in the Windows Control Panel) for machines running RayPlan should be set to "High Performance" at all times. This ensures that maximum processor performance is available for processor operations.

System monitoring

It is highly recommended to monitor the service operation, disk space availability, processor usage and additional key performance indicators on all RayPlan machines (Database and clients) with alerting of relevant resources where acceptable levels are breached. When parameters are outside acceptable levels, remedial actions should be assessed/performed as required. Services critical to the function of SQL and DICOM Storage SCP should be considered high priority for monitoring.

Whenever an issue is detected, there should always be a follow-up to identify the root cause following the resolution of the issue. Actions to address the root cause should be considered for implementation.

System Environment Acceptance Test

Whenever the RayPlan software or hardware is updated or in any way altered, the RayPlan System Environment Acceptance Test shall be performed (*RSL-D-RP-10A-SEAT, RayPlan 10A System Environment Acceptance Test Protocol*).

Disaster Recovery (DR) strategies

It is highly recommended that DR options are assessed, considered and where applicable implemented for RayPlan environments. These should be in line with the organization's policy for DR.

5.2 RAYPLAN STORAGE SCP SERVER

Where DICOM data is to be hosted by RaySearch software, the program can be run on the SQL server. The primary use case for the Storage SCP service is where a DICOM entity needs to "push" data to RayPlan. This is achieved by "pushing" to the service and (when configured in RayPlan) the data is available to the RayPlan users. Care should be taken regarding the port to be used by the service if any other software is installed on the machine.

5.3 THIRD PARTY SOFTWARE

Anti-Virus

Any Anti-Virus software installed on systems with RayPlan should have exemptions included for the RayPlan software. It is highly recommended that exceptions are applied on the RaySearch Laboratories folder within %programfiles% as well as all subfolders.

Windows updates

Windows updates should be assessed and applied in line with *RSL-D-RP-10A-IFU, RayPlan 10A Instructions for use*.

PDF Reader

A PDF reader software should be installed on every machine with RayPlan to enable users to read RayPlan documentation. An example would be Adobe Reader.

5.4 ACTIVE DIRECTORY (AD)

AD requirements and recommendations

Active Directory is required for all RayPlan installations with the exception of stand-alone machines. Where Active Directory is to be installed as part of an overall RayPlan installation (i.e., it is not present already) it should be noted that the installation of SQL Server and Active Directory on the same machine is not recommended by Microsoft or RaySearch. It is required to use a separate machine for Active Directory and SQL Server.

AD maintenance

It is highly recommended to test the authentication mechanisms in RayPlan following any Active Directory maintenance.

AD groups

The following AD groups are used for RayPlan operations that require additional authentication:

- **RayStation-Administration**¹ (access to Clinic Settings, ability to delete patients and other administrative tasks)
- **RayStation-BeamCommissioning**² (ability to commission machines within RayPhysics)
- **RayStation-PlanApproval**³ (ability to approve plans)

It is also recommended to create a RayPlan-Users group which can be used to provide basic access to RayPlan. The groups listed above can be nested to facilitate maintenance. Relevant subgroups can also be added for specific user profiles.

1 RayPlan-Administration or equivalent group name can be used with a group mapping in Clinic Settings for consistency if required.

2 RayPlan-BeamCommissioning or equivalent group name can be used with a group mapping in Clinic Settings for consistency if required.

3 RayPlan-PlanApproval or equivalent group name can be used with a group mapping in Clinic Settings for consistency if required.

AD Service accounts

Active Directory Service, Managed Service or Group Managed Service accounts can be considered for running SQL, Indexing, SCP or other services. The relevant file, folder, database or similar permissions should be granted to these service accounts.

AD Functional level

Windows Active Directory Functional level of Server 2012 R2 is required. RayPlan deployments not involving HIPAA related functionality can function within a Microsoft Active Directory Server 2008 R2 Functional level.

5.5 AD ORGANIZATIONAL UNIT

The creation of a specific Organisational Unit for RayPlan within an existing Domain should be considered as a best practice. This allows for AD policies to be tailored to the RayPlan environment.

5.6 SQL

SQL Backups

Details of how to set up scheduled SQL backups are available from RaySearch Laboratories. Existing SQL backup agents can be used by the clinic. All backup operations should be monitored and reviewed on a regular basis. Databases may grow at an accelerated rate where MS SQL internal processes are not triggered by MS SQL backups.

In RayPlan 10A the initial view of RayPlan (before opening a patient) provides alerts where full, partial or transaction log backups have not been successful for the system that RayPlan is connected to (Primary Patient DB only).

SQL Indexing

It is highly recommended that re-indexing of the RayPlan Databases occurs on a regular basis, e.g., following every full database backup.

Note: *Failure to re-index databases can lead to data corruption.*

SQL Access

When creating databases, all domain users have theoretical database access. This should be reviewed to ensure the correct user sets have access to the databases. When upgrading or creating additional databases access permissions should also be reviewed.

5.7 LICENSING OF THIRD PARTY SOFTWARE

Obtaining licenses for the Operating System and other supporting software such as SQL related to RayPlan is the responsibility of the customer clinic. The clinic should confirm that they have the correct numbers of applicable licenses with their local license vendor on a regular basis.

Licenses include (but are not limited to):

- Microsoft Server and Client Operating System Licenses – including annual fees under certain OS licenses

- Microsoft Server Client Access Licenses
- Microsoft SQL Server Client Access Licenses

5.8 RAM

RAM DIMMs should be installed to provide optimal RAM usage, i.e., units of 6 for Intel Xeon Scalable (48/96 GB) or units of 4 for Intel Xeon E5 (32/64 GB).

6 MICROSOFT REFERENCE INFORMATION

If you are not familiar with all Microsoft technologies, significant information is available from <https://mva.microsoft.com>.

Information related to RayPlan:

Active Directory: <https://mva.microsoft.com/en-US/training-courses/understanding-active-directory-8233>

SQL Server: <https://mva.microsoft.com/en-US/training-courses/sql-database-fundamentals-16944>



CONTACT INFORMATION



RaySearch Laboratories AB (publ)
Sveavägen 44
SE-111 34 Stockholm
Sweden

Contact details head office

P.O. Box 3297
SE-103 65 Stockholm, Sweden
Phone: +46 8 510 530 00
Fax: +46 8 510 530 30
info@raysearchlabs.com
www.raysearchlabs.com

RaySearch Americas

Phone: +1 877 778 3849

RaySearch France

Phone: +33 975 433 632

RaySearch Korea

Phone: +82 10 2230 2046

RaySearch Belgium

Phone: +32 2 213 83 65

RaySearch Germany

Phone: +49 30 89 36 06 90

RaySearch Singapore

Phone: +65 81 28 59 80

RaySearch China

Phone: +86 137 0111 5932

RaySearch Japan

Phone: +81 3 4405 6902

RaySearch UK

Phone: +44 7508 426 563